

CHEMISTRY LAB PREPARATION LIST FOR JC2 H2 CHEMISTRY

2018 PRELIM EXAM

APPARATUS NEEDED (PER STUDENT)

S/N	Apparatus	Quantity	Location	Remarks
1	10 cm ³ measuring cylinder	2	Student's bench	Reuse
2	50 cm ³ measuring cylinder	1	Student's bench	Reuse
3	25 cm ³ measuring cylinder	1	Student's bench	Reuse
4	100 cm ³ measuring cylinder	1	Student's bench	Reuse
5	stopwatch	1	Student's bench	Reuse
6	Thermometer (-10°C to +110°C)	1	Student's bench	Reuse
7	Glass rod	1	Student's bench	Reuse
8	Wash bottle of distilled water	1	Student's bench	Reuse
9	Pen suitable for labelling glassware	1	Student's bench	Reuse
10	Heat proof mat	1	Student's bench	Reuse
11	Filter Funnel	1	Student's bench	Reuse
12	Test tube rack	1	Student's bench	Reuse
13	Test tube holder	1	Student's bench	Reuse
14	Bunsen burner	1	Student's bench	Reuse
15	Tripod with gauze	1	Student's bench	Reuse
16	Lighter	1	Student's bench	Reuse
17	Spatula	1	Student's bench	Reuse
18	100 cm ³ conical flask	1	Student's bench	Reuse
19	Stand and bigger clamp	1	Student's bench	Reuse
20	plastic tub	1	Student's bench	trough for gas collection (minimum capacity 1 dm ³)
21	Paper towels		Student's bench	Must change to new one per shift
22	Test tubes	6	Student's bench	Must change to new one per shift
23	Boiling tubes	2	Student's bench	Must change to new one per shift
24	Filter Paper	2-3	Student's bench	Must change to new one per shift
25	Dropping Pipette	5	Student's bench	Must change to new one per shift

* access to a balance weighing

CHEMICALS NEEDED (PER STUDENT)

S/N	Chemical	Concentration	Quantity	Location	Remarks
1	Hydrochloric acid (FA 1)	2.0 mol dm ⁻³	70 cm ³		Dilute 170 cm ³ of concentrated (35-37%; approximately 11 mol dm ⁻³) hydrochloric acid to 1 dm ³
2	Magnesium (ribbon) (FA 2)	about 4 cm length	3 pieces		4 cm Magnesium ribbon supplied in a stoppered container labelled FA 2 .
3	Sodium thiosulfate (FA 3)	0.10 mol dm ⁻³	200 cm ³		Dissolve 24.82g of Na ₂ S ₂ O ₃ .5H ₂ O in each dm ³ of solution.
4	Hydrochloric acid (FA 4)	1.0 mol dm ⁻³	200 cm ³		Dilute 500 cm ³ of 2.0 mol dm ⁻³ HC/ to 1 dm ³ with distilled water.
5	manganese(IV) oxide, MnO ₂ (FA 5)	-	About 6 g		In a stoppered weighing bottle labelled FA 5;
6	sulfuric acid (FA 6)	0.5 mol dm ⁻³	40 cm ³		Labelled FA 6 Cautiously pour 27.5 cm ³ of concentrated (98%) sulfuric acid into 500 cm ³ of distilled water with continuous stirring. Make the solution up to 1 dm ³ with distilled water. Care: concentrated H ₂ SO ₄ is very corrosive [MH]
7	Solid sodium ethanedioate (FA 7)	-	about 2 g		In a stoppered bottle labelled FA 7
8	manganese(II) sulfate, MnSO ₄ FA 8)	About 0.5 mol dm ⁻³	3 cm ³		dissolving, in deionised water, about 84.5 g of MnSO ₄ .H ₂ O in each dm ³ of solution. This solution must be stored in a a stoppered vial labelled FA 8
9	cyclohexene, C ₆ H ₁₀ , (FA 9)	-	about 1 cm ³		in a stoppered vial labelled FA 9;
10	bromine water, Br ₂ (aq) (FA 10)	-	about 2 cm ³		in a stoppered vial labelled FA 10.

STANDARD BENCH REAGENTS

Hazard	Label	Identity	Notes (hazards given in this column are for the raw materials)
[MH]	Dilute hydrochloric acid	2.0 mol dm ⁻³ HCl	Dilute 170 cm ³ of concentrated (35 – 37%; approximately 11 mol dm ⁻³) hydrochloric acid to 1 dm ³
[C]	Dilute nitric acid	2.0 mol dm ⁻³ HNO ₃	Dilute 128 cm ³ of concentrated (70% w/v) nitric acid to 1 dm ³
[MH]	Dilute sulfuric acid	1.0 mol dm ⁻³ H ₂ SO ₄	Cautiously pour 55 cm ³ of concentrated (98%) sulfuric acid into 500 cm ³ of distilled water with continuous stirring. Make the solution up to 1 dm ³ with distilled water. Care: concentrated H ₂ SO ₄ is very corrosive
[C][MH][N]	Aqueous ammonia	2.0 mol dm ⁻³ NH ₃	Dilute 112 cm ³ of concentrated (35%) ammonia to 1 dm ³
[C]	Aqueous sodium hydroxide	2.0 mol dm ⁻³ NaOH	Dissolve 80.0 g of NaOH in each dm ³ of solution. Care: the process of solution is exothermic and any concentrated solution is very corrosive.
[MH]	Aqueous barium nitrate	0.1 mol dm ⁻³ barium nitrate	Dissolve 26.1 g of Ba(NO ₃) ₂ in each dm ³ of solution.
[N]	Aqueous silver nitrate	0.05 mol dm ⁻³ silver nitrate	Dissolve 8.5g of AgNO ₃ in each dm ³ of solution.
[MH]	Limewater	Saturated aqueous calcium hydroxide, Ca(OH) ₂	Prepare fresh limewater by leaving distilled water to stand over solid calcium hydroxide for several days, shaking occasionally. Decant or filter the solution.
	Aqueous potassium manganate (VII)	0.02 mol dm ⁻³ potassium manganate (VII)	Dissolve 3.16g of KMnO ₄ in each dm ³ of solution.
	Aqueous Iodine		Dissolve the 2.5g of solid in water containing sufficient potassium iodide to dissolve the iodine, then make up the solution to 1 dm ³ with distilled water.

The following materials and apparatus should be available **FOR EACH STUDENT**.

Red and blue litmus papers, plain filter strips for use with acidified manganate(VII), aluminium foil for testing nitrate/nitrite, wooden splints and the apparatus normally used in the Centre for use with limewater testing for carbon dioxide.

