

Confidential Instructions:

Candidates are advised to spend no more than:

- **60 minutes** on **Question 1**.
- **60 minutes** on **Question 2**.
- **30 minutes** on **Question 3**.

Preparation of materials**Question 1**

- i. 1.0 mol dm⁻³ sucrose solution may be prepared the day before the examination.
It should be kept in a covered container in a refrigerator.
0.5% methylene blue solution may be prepared the day before the examination.
It should be kept in a covered container.
The solutions must be at **room temperature** for the examination.
- ii. **S**, 1.0 mol dm⁻³ sucrose solution
This is prepared by sprinkling 68.4 g of sucrose, a little at a time, onto the surface of 80 cm³ of distilled water, stirring continuously as you sprinkle. Make up to 200 cm³ with distilled water.
- iii. **P**, at least 7 pieces of peeled potato wrapped in a damp paper towel in a covered dish, labelled **P**.
You may use any variety of the white (or Irish) potato, *Solanum tuberosum*.
Cut each piece of potato with a cross-sectional area of 1.5 cm × 1 cm.
Each candidate should be provided with a mixture of different lengths, varying from 4.5 cm to 6 cm.
The potato pieces for each candidate should be prepared on the day of the examination.
- iv. **M**, 0.5% methylene blue solution
This is prepared by putting 0.5 g of methylene blue into 80 cm³ of distilled water and stirring continuously. Make up to 100 cm³ with distilled water.

Question 2

- i. 15 cm³ of 0.1% 2,6-dichlorophenol indophenol (DCPIP) solution in a corked, labelled container (labelled as **DCPIP**). Dissolve powder in distilled water.
- ii. 5 cm³ of a standard solution of ascorbic acid (labelled as **vitamin C solution**) prepared as follows:
Dissolve 400 mg of ascorbic acid powder in 100 cm³ distilled water. Add 15 drops of BDH Universal Indicator solution and adjust the pH to between 7 and 8 by adding 5% sodium hydroxide solution drop by drop. This is **4.0mgcm⁻³** ascorbic acid solution.
The solution should be dispensed to students in a corked, **labelled** tube from which 1cm³ volumes can be withdrawn.

- iii. 5 cm³ each of fresh lemon juice, orange juice, grapefruit juice. Dispense to students in corked tubes **labelled** lemon juice, orange juice and grapefruit juice.

Apparatus List

Each candidate will require:

materials and apparatus for each candidate	quantity
1.0 mol dm ⁻³ sucrose solution in a beaker or container, labelled S , provided at room temperature (see Preparation of materials)	at least 200 cm ³
Pieces of potato in a beaker or container, labelled P (see Preparation of materials)	at least 7 pieces
Distilled water in a beaker or container, labelled W , provided at room temperature	at least 200 cm ³
Methylene blue solution, labelled M , provided at room temperature (see Preparation of materials)	at least 15 cm ³
10 cm ³ or 20 cm ³ syringes, with the means to wash them out	2
2 cm ³ or 3 cm ³ syringe, with the means to wash it out	1
Pipettes, plastic or glass with a teat	2
Beakers or containers (capacity 75 cm ³ to 100 cm ³)	5
Test-tubes – large (to hold more than 25 cm ³ but no more than 50 cm ³)	5
Test-tube rack to hold 5 large test-tubes	1
Test-tubes – small (capacity 20 cm ³ to 30 cm ³)	5
Test-tube rack to hold 5 small test-tubes	1
Glass rod	1
Ruler, marked in mm	1
Scalpel or sharp blade	1
White tile or surface for cutting	1
Container with tap water (capacity approximately 200 cm ³), labelled For washing	1
Container (capacity approximately 200 cm ³), labelled For waste	1
Paper towels	8
Glass marker pen	1
Stop-clock or timer showing seconds	1

Question 2

1	DCPIP	11	4 teat pipettes (droppers)
2	Vitamin C solution	12	4 50-ml beakers
3	Lemon juice	13	1 stopwatch
4	Orange juice	14	1 wash bottle labelled as distilled water
5	Grapefruit juice	15	1 permanent marker
6	6 Test tubes	16	1 white tile
7	1 test tube rack		
8	3 2-ml syringes		
9	2 5-ml syringes		
10	1 Glass rod		

Question 3

- 1 Microscope with an eyepiece graticule fitted into the eyepiece lens
For each candidate:
 - the microscope must be set up on low power
 - the slide must not be left on the stage of the microscope.
- 2 Slide **M1** – *Ranunculus t.s.* root metaxylem