



NATIONAL JUNIOR COLLEGE
SH2 Year-End Practical Examination
Higher 2

CANDIDATE
NAME

SUBJECT
CLASS

REGISTRATION
NUMBER

CHEMISTRY

9729/04

Paper 4 Practical

16 August 2022

Candidates answer on the Question paper

2 hours 30 minutes

Additional Materials: As listed in the Confidential Instructions

READ THESE INSTRUCTIONS FIRST

Write your identification number and name.

Give details of the practical shift and laboratory where appropriate, in the boxes provided.

Write in blue or black pen.

You may use an HB pencil for any diagrams or graphs.

Do not use staples, paper clips, glue or correction fluid.

Answer **all** questions in the spaces provided on the Question Paper.

The use of an approved scientific calculator is expected, where appropriate. You may lose marks if you do not show your working or if you do not use appropriate units.

Qualitative Analysis Notes are printed on pages 22 and 23.

At the end of the examination, fasten all your work securely together.

The number of marks is given in brackets [] at the end of each question or part question.

Shift
Laboratory

For Examiner's use	
1	/ 9
2	/ 21
3	/ 17
4	/ 8
Total	/ 55

This document consists of **23** printed pages including this cover page.

Answer **all** the questions in the spaces provided.


1 Determination of the enthalpy change of a reaction, ΔH ,

FA 1 is 1.00 mol dm⁻³ sodium hydrogen carbonate, NaHCO₃ (**to be used for Q3 as well**)

FA 2 is 2.00 mol dm⁻³ sodium hydroxide, NaOH

FA 3 is 2.00 mol dm⁻³ sulfuric acid, H₂SO₄


An acid-base neutralisation reaction involves reacting the two solutions, to produce water molecules. The equation for this neutralisation reaction is given below.



You will follow the instructions to perform two experiments, **Experiment A** and **Experiment B**. Record your results in **Tables 1.1** and **1.2**.

Experiment A

Reaction between **FA 1**, NaHCO₃, and **FA 2**, NaOH.



The molar enthalpy change for **reaction 1**, $\Delta H_{\text{reaction1}}$, is the enthalpy change when 1.00 mol of NaHCO₃ reacts completely with NaOH.

- Using a 50 cm³ measuring cylinder, transfer 30.0 cm³ of **FA 1** into a Styrofoam cup. Place this cup inside a second Styrofoam cup, which is placed in a 250 cm³ glass beaker. Place the lid on the cup.
- Stir and measure the temperature of this **FA 1**, T_{FA1} .
- Using another 50 cm³ measuring cylinder, measure 20.0 cm³ of **FA 2**.
- Stir and measure the temperature of this **FA 2**, T_{FA2} .
- Add **FA 2** from the measuring cylinder to the **FA 1** in the Styrofoam cup. Immediately replace the lid.
- Using the thermometer, stir the mixture continuously until a maximum temperature is reached. Read and record this temperature T_{max} .
- Calculate the weighted average initial temperature, T_{average} , of **FA 1** and **FA 2** using the formula given below:


$$T_{\text{average}} = \frac{(V_{\text{FA1}} \times T_{\text{FA1}}) + (V_{\text{FA2}} \times T_{\text{FA2}})}{(V_{\text{FA1}} + V_{\text{FA2}})}$$

	Experiment A
$T_{\text{FA1}} / ^\circ\text{C}$	
$T_{\text{FA2}} / ^\circ\text{C}$	
$T_{\text{average}} / ^\circ\text{C}$	
$T_{\text{max}} / ^\circ\text{C}$	
$\Delta T_{\text{max}} / ^\circ\text{C}$	

Table 1.1

Experiment B

Reaction between **FA 1**, NaHCO_3 , and **FA 3**, H_2SO_4 .



The molar enthalpy change for **reaction 2**, $\Delta H_{\text{reaction 2}}$, is the enthalpy change when 1.00 mol of NaHCO_3 reacts completely with H_2SO_4 .

- 1 Using a measuring cylinder, transfer 30.0 cm³ of **FA 1** into a Styrofoam cup. Place this cup inside a second Styrofoam cup, which is placed in a 250 cm³ glass beaker. Place the lid on the cup.
- 2 Stir and measure the temperature of this **FA 1**, T_{FA1} .
- 3 Using another measuring cylinder, measure 20.0 cm³ of **FA 3**.
- 4 Stir and measure the temperature of this **FA 3**, T_{FA3} .
- 5 Add **slowly**, the **FA 3** from the measuring cylinder to the **FA 1** in the Styrofoam cup. Immediately replace the lid.
- 6 Using the thermometer, stir the mixture continuously until a minimum temperature is reached. Read and record this temperature T_{min} .
- 7 Calculate the weighted average initial temperature, T_{average} , of **FA 1** and **FA 3** using the formula given below:

$$T_{\text{average}} = \frac{(V_{\text{FA1}} \times T_{\text{FA1}}) + (V_{\text{FA3}} \times T_{\text{FA3}})}{(V_{\text{FA1}} + V_{\text{FA3}})}$$

	Experiment B
$T_{\text{FA1}} / ^\circ\text{C}$	
$T_{\text{FA3}} / ^\circ\text{C}$	
$T_{\text{average}} / ^\circ\text{C}$	
$T_{\text{min}} / ^\circ\text{C}$	
$\Delta T_{\text{max}} / ^\circ\text{C}$	

Table 1.2

[2]

- (a) For the purpose of calculations, you should assume that the mixture has a density of 1.00 g cm^{-3} and specific heat capacity, c , of $4.18 \text{ J g}^{-1} \text{ K}^{-1}$.
- (i) Use your results from **Table 1.1** to calculate a value for the molar enthalpy change for **reaction 1**, $\Delta H_{\text{reaction1}}$.

$$\Delta H_{\text{reaction1}} = \dots\dots\dots [2]$$


- (ii) Use your results from **Table 1.2** to calculate a value for the molar enthalpy change for **reaction 2**, $\Delta H_{\text{reaction2}}$.

$$\Delta H_{\text{reaction2}} = \dots\dots\dots [2]$$

- (b) Ionic equations for neutralisation, **reaction 1**, and **reaction 2** are shown below.

$\text{H}^+(\text{aq}) + \text{OH}^-(\text{aq}) \rightarrow \text{H}_2\text{O}(\text{l})$	$\Delta H_{\text{neu}} = -57.1 \text{ kJ mol}^{-1}$
$\text{HCO}_3^-(\text{aq}) + \text{OH}^-(\text{aq}) \rightarrow \text{CO}_3^{2-}(\text{aq}) + \text{H}_2\text{O}(\text{l})$	$\Delta H_{\text{reaction1}}$
$\text{HCO}_3^-(\text{aq}) + \text{H}^+(\text{aq}) \rightarrow \text{H}_2\text{O}(\text{l}) + \text{CO}_2(\text{g})$	$\Delta H_{\text{reaction2}}$

Carbon dioxide reacts with solutions of carbonate ions according to the following equation.



Using your calculated answers in (a), together with the given value of enthalpy change of neutralisation, ΔH_{neu} , construct an energy cycle to determine a value for the enthalpy change for this reaction, $\Delta H_{\text{reaction3}}$.

$\Delta H_{\text{reaction3}} = \dots\dots\dots$
[3]

[Total: 9]

2 Qualitative analysis of some organic and inorganic compounds

You are provided with liquids **FA 4**, **FA 5**, **FA 6**, **FA 7**, and **FA 8**.

You are to perform the tests described in **Tables 2.1** and **2.3** and record your observations.

At each stage of any test, you are to record details of the following:

- colour changes seen
- the formation of any precipitate and its solubility in an excess of the reagent added
- the formation of any gas and its identification by a suitable test
- if there is no observable change, write **no observable change**

You should indicate clearly at which stage in a test a change occurs, recording your observations alongside the relevant tests.

In all tests, the reagents should be added gradually until no further change is observed unless you are instructed otherwise. If any solution is warmed, a **boiling tube** must be used. Rinse and reuse test-tubes where possible.

No additional or confirmatory tests for ions present should be attempted.

(a) **FA 4**, **FA 5**, **FA 6** and **FA 7** are organic compounds with the molecular formula shown.

- **FA 4**: $\text{C}_3\text{H}_6\text{O}$
- **FA 5**: $\text{C}_3\text{H}_6\text{O}$
- **FA 6**: $\text{C}_3\text{H}_8\text{O}$
- **FA 7**: $\text{C}_3\text{H}_6\text{O}_2$

You will perform some of the tests described in **Table 2.1**.

Using the observations in **Table 2.1** and the given molecular formula, you will then deduce the identities of **FA 4**, **FA 5**, **FA 6** and **FA 7**.

Perform the tests given in **Table 2.1**. Some of the observations have been completed for you. There is no need to carry out those tests. Record your observations in **Table 2.1**.

Safety: Organic compounds are flammable. Transfer your organic waste into the waste bottle for disposal after the end of the assessment.

Test	Procedure	Observations with FA 4, C ₃ H ₆ O	Observations with FA 5, C ₃ H ₆ O	Observations with FA 6, C ₃ H ₈ O	Observations with FA 7, C ₃ H ₆ O ₂
1.	[2,4-DNPH] To 1 cm depth of unknown liquid, add 1cm depth of 2,4-dinitrophenylhydrazine solution.		Orange precipitate observed		
2.	To 1 cm depth of unknown liquid, add aqueous sodium carbonate. Test for and identify any gas evolved.	No observable change	No observable change	No observable change	Effervescence produced. White precipitate observed when gas is bubbled into limewater
3.	[Iodoform Test] To 1 cm depth of unknown liquid, add 8 drops of aqueous sodium hydroxide followed by iodine solution, dropwise, until a permanent orange/red colour is present. Warm the mixture in a beaker of hot water for two minutes.	No observable change	Yellow precipitate observed		No observable change

4.	[Tollens' reagent] To 1 cm depth of aqueous silver nitrate, slowly add 1 cm depth of aqueous sodium hydroxide. Add aqueous ammonia slowly, with shaking, until the precipitate formed dissolves. You may use a clean glass rod to stir the mixture and help dissolve the precipitate. Add 1 cm depth of unknown liquid to this mixture, shake the tube and let it stand.				
----	--	--	--	--	--

Table 2.1

[4]

(b) Complete **Table 2.2** with the identities of **FA 4**, **FA 5**, **FA 6**, and **FA 7**.

Give evidence from the observations in **Table 2.1** to support your conclusions.

	Compound	Evidence
FA 4		
FA 5		
FA 6		
FA 7		

Table 2.2

[4]

FA 8 is a mixture of two salts and contains two cations.

You will perform a series of test-tube reactions and use the observations to help you identify the two cations.

(c)	Test	Observation
1.	Test the FA 8 solution using Universal indicator paper.	
2.	Prepare about 2 cm³ of dilute sulfuric acid in a clean test tube. Place about 2 cm³ of FA 8 into another clean test-tube. Carefully add aqueous sodium hydroxide, dropwise with shaking, until no further change is seen. Swirl and filter the mixture, collecting the filtrate in the test tube containing dilute sulfuric acid. Observe the residue and filtrate until no further change is seen.	
3.	To 1cm ³ of FA 8 , carefully add aqueous potassium manganate, dropwise with shaking, until no further change is seen.	

Table 2.3

(d) (i) Identify the cations in **FA 8**.

[5]

Cation 1:

Cation 2:

[2]

(ii) With the aid of an equation, explain the observation for **Test 1** in **Table 2.3**.

.....

.....

.....

[2]

- (iii) Identify the metal-containing complex formed when excess sodium hydroxide is added to **FA 8**.

Write equations to illustrate the formation of this complex.

.....

.....

.....

[2]

- (iv) **FA 8** contains either the sulfate ion or a halide anion.

Describe two tests, using only the bench reagents provided, which will allow you to identify the anion present.

In each case, state how you will decide if the test result is positive.

[You **DO NOT** need to carry out these tests.]

.....

.....

.....

.....

.....

.....

.....

.....

[2]

[Total: 21]

3 To determine the order of reaction with respect to the concentration of iodine in the iodination of propanone reaction

FA 12 is 1.00 mol dm⁻³ propanone, CH₃COCH₃

FA 13 is 1.00 mol dm⁻³ sulfuric acid, H₂SO₄


FA 14 is an aqueous solution of iodine, I₂

FA 15 is 0.0150 mol dm⁻³ sodium thiosulfate, Na₂S₂O₃

FA 1 is 1.00 mol dm⁻³ sodium hydrogen carbonate, NaHCO₃ (**same solution used in Q1**)

You are also provided with a starch indicator.

The iodination of propanone, to form iodopropanone, proceeds as shown in the equation below.




This reaction is first order with respect to both CH₃COCH₃ and H⁺ ions.

You are to investigate the order of reaction with respect to I₂.

A reaction mixture containing **FA 12**, **FA 13**, and **FA 14** is first prepared. At different chosen times, aliquots (fixed volumes) of this reaction mixture is removed and quenched using excess **FA 1**.

It is necessary that you titrate each aliquot against **FA15** before transferring the next aliquot. I₂ and S₂O₃²⁻ react as shown in **Reaction 2**.



The required order of reaction can be obtained by the graphical analysis of your results.

The first aliquot should be removed approximately 4 minutes after the reagents were mixed. You will then remove four further aliquots, at time intervals of your choice, up to a maximum time of 30 minutes.

In an appropriate format in the space provided on **page 14**, prepare a table in which to record for each aliquot

- the time of transfer, *t*, in minutes and seconds,
- the decimal time, *t_d*, in minutes, to 0.1 min, for example, if *t* = 4 min 33 s then *t_d* = 4 min + 33/60 = 4.6 min,
- the burette readings and the volume of **FA 15** added.

Safety: Propanone is flammable. Transfer your titrated solutions into the waste bottle for later disposal.

Keep the conical flask (**reaction mixture**) stoppered except when removing aliquots.

Experiment

NOTE READ THE FULL PROCEDURE BEFORE STARTING YOUR EXPERIMENT

1. Label each of the boiling tubes **1** to **5**.
2. Add approximately 10 cm³ of **FA 1** to each of these boiling tubes.
3. Fill a burette with **FA 15**

Preparing the reaction mixture

4. Using a 50 cm³ measuring cylinder, transfer 25.0 cm³ of **FA 12** into the 100 cm³ beaker.
5. Using the same 50 cm³ measuring cylinder, transfer 25.0 cm³ of **FA 13** into the same 100 cm³ beaker.
6. Using another 50 cm³ measuring cylinder, transfer 50.0 cm³ of **FA 14** into the 250 cm³ conical flask, labelled **reaction mixture**.
7. Pour the contents of the 100 cm³ beaker into this 250 cm³ conical flask. Start the stopwatch, **insert the stopper** and swirl the mixture thoroughly. Once you have started the stopwatch, allow it to continue running for the duration of the experiment. You must not stop the stopwatch until you have collected all of your aliquots.

Removing aliquots of reaction mixture and titration

8. At approximately 4 minutes, using a 10.0 cm³ pipette, remove a 10.0 cm³ aliquot of the reaction mixture.
9. **Immediately** transfer this aliquot into the boiling tube labelled **1** and swirl the mixture. Read and record the transfer time in minutes and seconds, to the nearest second, when half of the reaction mixture has emptied from the pipette. Replace the stopper in the flask.
10. Pour all the contents of **boiling tube 1** into a clean 250 cm³ conical flask. Wash out the boiling tube and add the washings to the conical flask.
11. Titrate the iodine in this solution with **FA 15**. Add about 1 cm³ of starch indicator when the colour of the solution turns pale yellow. The solution will turn blue-black. The end-point is reached when the dark blue-black colour just disappears. Record your results.
12. Empty the contents of this conical flask into the waste bottle. Wash this conical flask thoroughly with water.
13. At approximately 8 minutes, repeat points **8** to **12**. Transfer this aliquot into the boiling tube labelled **2**.
14. Repeat point **8** to **12** for the remaining boiling tubes at about four-minute intervals.

(a) Results

[4]

- (b) (i) On **Fig 3.1**, plot a graph of **volume of sodium thiosulfate, FA 15**, on the y-axis, against decimal time, t_d , on the x-axis. Start the x-axis at $t_d = 0$. You should choose a scale which will allow you to extrapolate your graph back to $t_d = 0$.

Draw the most appropriate best-fit line taking into account all of your plotted points.




Fig 3.1

[3]

- (ii) Deduce the order of reaction with respect to the I_2 in reaction 1. Explain your answer.

order.....

explanation.....

.....

.....

[1]

- (c) (i) Write the rate equation for the iodination of propanone.

..... [1]

- (ii) Calculate the gradient of the line you have drawn in **Fig 3.1**, showing clearly how you did this.

gradient = $\text{cm}^3 \text{min}^{-1}$

[1]

- (iii) Use your answer from (c)(ii) to determine the rate of change of amount of $S_2O_3^{2-}$ ions required in mol min^{-1} .

rate of change of amount of $S_2O_3^{2-}$ ions required = mol min^{-1}

[1]

- (iv) Hence, deduce the rate of disappearance of I_2 in mol min^{-1} .

rate of disappearance of I_2 = mol min^{-1}

[1]

- (v) Use your answer from (c)(iv) to calculate the rate of change of $[I_2]$ in the reaction mixture.

rate of change of $[I_2]$ in the reaction mixture = mol dm⁻³ min⁻¹
[1]

- (vi) Hence, calculate the value of the rate constant for this reaction, giving its units.

rate constant =
[2]

- (d) **Step 7** requires you to mix each aliquot immediately with an excess of sodium hydrogencarbonate solution, **FA 1**. Suggest a clear explanation for this requirement.

.....
.....
[1]


- (e) Explain why the concentration of iodine in **FA 14** used is very much lower than the concentrations of propanone in **FA 12** and of H^+ ions in **FA 13**.

.....
.....
.....
.....
[1]

[Total: 17]

4 Planning – Extension to Q3

You are to plan a series of experiments to verify that the order with respect to propanone in the iodination of propanone is 1.



This investigation can be carried out by monitoring the time taken for the colour of a limited amount of iodine to be discharged when reacting with varying concentrations of propanone.

You may assume that it takes approximately 1 minute for the iodine colour to be discharged when a mixture of 20 cm³ of 2.00 mol dm⁻³ of propanone, CH₃COCH₃ and 10 cm³ of 1.00 mol dm⁻³ sulfuric acid, H₂SO₄ is reacted with 10 cm³ of 1.00 mol dm⁻³ iodine, I₂.

You may assume that you are provided with:

- **Only 60 cm³** of 2.00 mol dm⁻³ propanone, CH₃COCH₃
- 1.00 mol dm⁻³ sulfuric acid, H₂SO₄
- 1.00 mol dm⁻³ solution of iodine, I₂
- 100 cm³ volumetric flasks
- starch indicator
- stopwatch
- the equipment normally found in a school or college laboratory.

(a) Serial dilution

Suggest how propanone of concentrations 1.00 mol dm⁻³, 0.50 mol dm⁻³ 0.25 mol dm⁻³ and 0.125 mol dm⁻³ can be prepared from the 2.00 mol dm⁻³ propanone provided.

Clearly outline the procedure for the dilution.

.....

.....

.....

.....

.....

.....

.....

.....

[2]

- (b) Plan a procedure to collect sufficient data to allow a graph to be plotted to verify the order of reaction with respect to propanone.

You are to use the 1.00 mol dm^{-3} , 0.50 mol dm^{-3} , 0.25 mol dm^{-3} and $0.125 \text{ mol dm}^{-3}$ propanone solutions you have prepared in (a) in your procedure.

In your plan you should include brief details of:

- the apparatus you would use,
- the quantity of the different solutions you would use,
- the procedure you would follow,
- how you would verify that the order with respect to propanone is 1.

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

- (c) A student suggested that preparing individual solutions of propanone of varying concentrations for use in this analysis is very time consuming and not necessary.

Suggest how the above procedure can be modified to investigate the order with respect to propanone by just using 2.00 mol dm^{-3} propanone.

You may find it helpful to use of a table to illustrate your answer.

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

[2]

[Total: 8]

Qualitative Analysis Notes

[ppt. = precipitate]

(a) Reactions of aqueous cations

cation	reaction with	
	NaOH(aq)	NH₃(aq)
aluminium, Al ³⁺ (aq)	white ppt. soluble in excess	white ppt. insoluble in excess
ammonium, NH ₄ ⁺ (aq)	ammonia produced on heating	–
barium, Ba ²⁺ (aq)	no ppt. (if reagents are pure)	no ppt.
calcium, Ca ²⁺ (aq)	white ppt. with high [Ca ²⁺ (aq)]	no ppt.
chromium(III), Cr ³⁺ (aq)	grey-green ppt. soluble in excess giving dark green solution	grey-green ppt. insoluble in excess
copper(II), Cu ²⁺ (aq)	pale blue ppt. insoluble in excess	blue ppt. soluble in excess giving dark blue solution
iron(II), Fe ²⁺ (aq)	green ppt., turning brown on contact with air insoluble in excess	green ppt., turning brown on contact with air insoluble in excess
iron(III), Fe ³⁺ (aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess
magnesium, Mg ²⁺ (aq)	white ppt. insoluble in excess	white ppt. insoluble in excess
manganese(II), Mn ²⁺ (aq)	off-white ppt., rapidly turning brown on contact with air insoluble in excess	off-white ppt., rapidly turning brown on contact with air insoluble in excess
zinc, Zn ²⁺ (aq)	white ppt. soluble in excess	white ppt. soluble in excess

(b) Reactions of anions

<i>anions</i>	<i>reaction</i>
carbonate, CO_3^{2-}	CO_2 liberated by dilute acids
chloride, Cl^- (aq)	gives white ppt. with Ag^+ (aq) (soluble in NH_3 (aq))
bromide, Br^- (aq)	gives pale cream ppt. with Ag^+ (aq) (partially soluble in NH_3 (aq))
iodide, I^- (aq)	gives yellow ppt. with Ag^+ (aq) (insoluble in NH_3 (aq))
nitrate, NO_3^- (aq)	NH_3 liberated on heating with OH^- (aq) and Al foil
nitrite, NO_2^- (aq)	NH_3 liberated on heating with OH^- (aq) and Al foil; NO liberated by dilute acids (colourless $\text{NO} \rightarrow$ (pale) brown NO_2 in air)
sulfate, SO_4^{2-} (aq)	gives white ppt. with Ba^{2+} (aq) (insoluble in excess dilute strong acids)
sulfite, SO_3^{2-} (aq)	SO_2 liberated by dilute acids; gives white ppt. with Ba^{2+} (aq) (soluble in dilute strong acids)

(c) Tests for gases

<i>gas</i>	<i>test and test result</i>
ammonia, NH_3	turns damp red litmus paper blue
carbon dioxide, CO_2	gives white ppt. with limewater (ppt. dissolves with excess CO_2)
chlorine, Cl_2	bleaches damp litmus paper
hydrogen, H_2	"pops" with a lighted splint
oxygen, O_2	relights a glowing splint
sulfur dioxide, SO_2	turns aqueous potassium manganate(VII) from purple to colourless

(d) Colour of halogens

<i>halogen</i>	<i>colour of element</i>	<i>colour in aqueous solution</i>	<i>colour in hexane</i>
chlorine, Cl_2	greenish yellow gas	pale yellow	pale yellow
bromine, Br_2	reddish brown gas / liquid	orange	orange-red
iodine, I_2	black solid / purple gas	brown	purple